

# CDC 3800 Computer Program to Calculate and Plot the Compressibility, Modulus of Compressibility, and Perpendicular Component of the Molecular Dipole Moment of Films Spread on Liquid Surfaces

LAURA A. HALPER

*Laboratory for Chemical Physics*

December 31, 1969



**NAVAL RESEARCH LABORATORY**  
**Washington, D.C.**

## ABSTRACT

A program of interest to surface chemists has been developed for use on NRL's Control Data Corporation 3800 computer. The program reads film pressure, surface potential, and molecular area data. Tabulated values are returned for the compressibility, modulus of compressibility, and perpendicular component of the molecular dipole moment, along with corresponding molecular areas. Library routines for the Calcomp Plotter are used to plot the compressibility vs area/molecule, modulus of compressibility vs area/molecule, and perpendicular component of the molecular dipole moment vs area/molecule.

## PROBLEM STATUS

This is an interim report; work is continuing on the problem.

## AUTHORIZATION

NRL Problem C02-10  
Project RR 001-01-43-4751

Manuscript submitted October 10, 1969.

CDC 3800 COMPUTER PROGRAM TO CALCULATE AND PLOT THE  
COMPRESSIBILITY, MODULUS OF COMPRESSIBILITY, AND  
PERPENDICULAR COMPONENT OF THE MOLECULAR  
DIPOLE MOMENT OF FILMS SPREAD ON LIQUID SURFACES

## INTRODUCTION

Much of the data collected in the study of films spread on liquid surfaces is in the form of either film pressure or surface potentials at various molecular areas. The interpretation of these data may indicate conformation and orientation of molecules in the liquid surface.

This program, titled "Films," was developed to avoid the tedious slope calculations necessary for obtaining compressibility curves from film pressure data. Because surface potential,  $\Delta V$ , data are often collected along with film pressure data, the program uses  $\Delta V$  data to calculate molecular dipole moments perpendicular to the liquid surface.

## GENERAL PROGRAM INFORMATION

Input is grouped either as set 1 data or set 2 data, depending on whether it is film pressure vs area/molecule (set 1) or surface potential vs area/molecule (set 2). For set 1 data, the program returns tabulations and plots of the compressibility and modulus of compressibility vs area/molecule. For set 2 data, the program returns perpendicular dipole vs area/molecule.

The program may be run with either or both sets present in the data deck. All set 1 data must be kept together, and all set 2 data must be kept together. Set 1 data need not precede set 2 data in the data deck. Input may be presented in either increasing or decreasing order by area, but successive data points must be presented successively in the data deck.

The program accepts up to and including 50 points for each film pressure run and 50 points for each surface potential run. The program accepts up to and including 99 runs for each set 1 and set 2.

Output quantities are calculated as follows:

$$\text{Compressibility}_{(i)} \text{ in cm/dyne} = \frac{1}{A_{(i)}} \cdot \frac{A_{(i-1)} - A_{(i)}}{F_{(i)} - F_{(i-1)}},$$

where  $(i)$  refers to the  $i$ th data point,  $F$  is the film pressure in dynes/cm, and  $A$  is the molecular area in square angstroms/molecule.

$$\text{Modulus}_{(i)} \text{ of compressibility in dynes/cm} = \frac{1}{\text{compressibility}_{(i)}}$$

and

$$\text{Perpendicular dipole}_{(i)} \text{ in Debye units} = \frac{\Delta V_{(i)} A_{(i)}}{4\pi} \times 3 \times 10^{-4},$$

where  $\Delta V$  is surface potential in millivolts.

#### INPUT LIMITS

Area/molecule: 0.00 to 999.99 square angstroms/molecule

Film pressure: 0.00 to 99.99 dynes/cm

Surface potential: 0.00 to 999.99 millivolts

#### OUTPUT LIMITS

Compressibility: no limits on magnitude but three significant digits

Modulus of compressibility: same as for compressibility

Perpendicular dipole: 0.00 to 99.99 Debye units

#### DEFINITIONS OF TERMS IN PROGRAM

ITSETS — the total number of data sets in the data deck. ITSETS will equal 1 if only one type of data (either set 1 or set 2) is present. ITSETS will equal 2 if both sets 1 and 2 are present in the data deck.

ITSET — designation of data type: ITSET = 1 indicates film pressure-vs-area data and ITSET = 2 indicates surface potential-vs-area data

ITRUNS — total number of experimental trials for either set 1 or set 2 data. For example, if one film-pressure-vs-area curve was made of five different compounds, ITRUNS for set 1 would equal 5. If five film pressure-vs-area curves were made of one compound, ITRUNS for set 1 would equal 5. If three film pressure-vs-area curves were made of one compound and one film pressure-vs-area curve was made of two additional compounds, ITRUNS for set 1 would equal 5.

$$0 \leq \text{ITRUNS} \leq 99$$

IRUN — ordinal number of an experimental trial in either set 1 or set 2.

$$0 \leq \text{IRUN} \leq 99$$

NUM — number of data points in some particular IRUN.

$$0 \leq \text{NUM} \leq 50$$

Information about the Calcomp Plotter subroutines used in this program can be found in the following notes and memoranda:

1. Memorandum 7810-5:ABB:rmd:pj, September 18, 1967, Plotter Subroutine Package.
2. Memorandum 4500-31:DEG:rmd, August 24, 1967, 3800 Utility Program.

3. Memorandum 7810-33:GHR:rmd, November 22, 1967, Plotter Subroutine Package.
4. Memorandum 7810-121:GR:vjs, August 14, 1968, Changes in Plotter Subroutine Package.
5. NRL Computer Note 6, October 17, 1968, Plot Time Message for the Calcomp Plotter.
6. NRL Computer Note 10, December 17, 1968, Initialization for the 3800 Calcomp Plotter Package.
7. NRL Computer Note 25, June 12, 1969, New SCALE Subroutine to go on the 3800 System Library.

## PROGRAM

```

PROGRAM FILMS
  DIMENSION AREA (50),PRESS(50),COMP(50),CMOD(50),DELV(50),PERMU(50)
1,DELA(50),DELP(50),ARRAY(254)
  CALL PLOTS (ARRAY,254,7)
  ICOUNT=0
  READ 100,ISETS
  READ 101,ISET,ITRUNS
8  GO TO (1,2) ISET
1  DO 4 I=1,ITRUNS
  READ 102,IRCN,NUM,N,A,M,E
  NNUM=NUM-1
  DO 5 II=1,N,M
  AREA(II)=0,
  PRESS(II)=0,
5  CONTINUE
  DO 6 II=1,N,M
  READ 103, AREA(II),PRESS(II)
6  CONTINUE
  DO 7 II=2,N,M
  DELA(II)=AREA(II-1)-AREA(II)
  DELP(II)=PRESS(II)-PRESS(II-1)
  COMP(II)=(1./AREA(II))*(DELA(II)/DELP(II))
  CMOD(II)=1./COMP(II)
7  CONTINUE
  PRINT 200,ISET,IRUN,N,A,M,E
  DO 9 J=1,NUM
  PRINT 201,PRESS(J),AREA(J)
9  CONTINUE
  PRINT 202
  DO 10 J=2,N,M
  PRINT 203,COMP(J),CMOD(J),AREA(J)
10 CONTINUE
  CALL SCALE(COMP(2),NNUM,8,,CMIN,DC,1,CTK)
  CALL AXIS(0,,0,,25HCOMPRESSIBILITY (CM/DYNE),25,8,,90,,CTK.CMIN,
1DC,5HE10,,3)
  CALL SCALE(AREA(2),NNUM,10,,AMIN,DA,1,ATK)
  CALL AXIS(0,,0,,28HAREA/MOLECULE (SQ.ANGSTROMS),-28,10,,0,,ATK,
1AMIN,DA,4HF6,2)
  CALL PLOT (0,,0,,2)
  CALL LINE (AREA(2),COMP(2),NNUM,1,2,0,1,1)
  CALL PLOT(1,,10,,3)
  CALL PLOT(1,,10,,2)
  CALL SYMBOL (1,,10,,0,21,21HCOMPRESSIBILITY CURVE,0,,21)
  CALL PLOT(13,,0,,3)
  CALL SCALE(CMOD(2),NNUM,8,,CMMIN,DCM,1,CMTK)
  CALL AXIS(0,,0,,17HMODULUS (DYNE/CM),17,8,,90,,CMTK,CMMIN,DCM,

```

```

15HE10,3)
CALL AXIS (0,,0,,28HAREA/MOLECULE (SQ,ANGSTROMS),-28,10,,0,,ATK,
1AMIN,DA,4HF6,2)
CALL PLOT(0,,0,,2)
CALL LINE (AREA(2),CMOD(2),NNU,1,1,0,1,1)
CALL PLOT(1,,10,,3)
CALL PLOT(1,,10,,2)
CALL SYMBOL (1,,10,,0,21,26HMODULUS OF COMPRESSIBILITY,0,,26)
CALL PLOT (13,,0,,3)
4   CONTINUE
ICOUNT=ICOUNT+1
20  IF (ITSETS, EQ, 1) GO TO 500
    IF (ICOUNT, EQ, 2) GO TO 500
    READ 101, ISET, ITRUNS
    GO TO 8
2   DO 15 I=1, ITRUNS
    READ 102, IRUN, NUM, N, A, M, E
    DO 16 II=1, NUM
    AREA (II)=0
    DELV(II)=0
10  CONTINUE
    DO 13 II=1, NUM
    READ 103, AREA(II), DELV(II)
    PERMU(II)=(DELV(II)*AREA(II))/(4,*3,14159)*3,* (10,**(4))
15  CONTINUE
    PRINT 204, ISET, IRUN, N, A, M, E
    DO 11 II=1, NUM
    PRINT 201, DELV(II), AREA(II)
11  CONTINUE
    PRINT 205
    DO 12 II=1, NUM
    PRINT 206, PERMU(II), AREA(II)
12  CONTINUE
    CALL SCALE (PERMU, NUM, 8,, PMIN, DP, 1, PTK)
    CALL AXIS (0,,0,,15HDIPOLE (DEBYES),15,8,,90,,PTK, PMIN, DP, 4HF4, 2)
    CALL SCALE (AREA, NUM, 10,, AMIN, DA, 1, ATK)
    CALL AXIS (0,,0,,28HAREA/MOLECULE (SQ,ANGSTROMS),-28,10,,0,,ATK,
1AMIN,DA,4HF6,2)
    CALL PLOT (0,,0,,2)
    CALL LINE (AREA, PERMU, NUM, 1, 0, 0, 1, 1)
    CALL PLOT (1,,10,,3)
    CALL PLOT (1,,10,,2)
    CALL SYMBOL (1,,10,,0,21,20HPERPENDICULAR DIPOLE,0,,20)
    CALL PLOT (13,,0,,3)
15  CONTINUE
    ICOUNT=ICOUNT+1
    GO TO 20
100  FORMAT (I1)
101  FORMAT (I1,5X,I2)
102  FORMAT (I2,5X,I2,5X,A8,A8,A8,A8)
103  FORMAT (F6,2,5X,F6,2)
200  FORMAT(///,27X,20HORIGINAL DATA: SET=,I1,2X,4HRUN=,I2,2X,4A8/
138X,23HFILM PRESSURE(DYNES/CM),5X,28HAREA/MOLECULE(SQ, ANGSTROMS)/
2)
201  FORMAT (47X,F6,2,25X,F6,2)
202  FORMAT (///,58X,15HCALCULATED DATA/17X, 24HCOMPRESSIBILITY(CM/DYNE
1),5X,35HMODULUS OF COMPRESSIBILITY(DYNE/CM),5X,28HAREA/MOLECULE(SQ
2, ANGSTROMS)/)
203  FORMAT (24X,E10,3,25X,E10,3,30X,F6,2)
204  FORMAT(///,28X,20HORIGINAL DATA: SET=,I1,2X,4HRUN=,I2,2X,4A8/
137X,25HDELTA VOLTAGE(MILLIVOLTS),5X,28HAREA/MOLECULE(SQ, ANGSTROMS
2)/)
205  FORMAT (///,58X,15HCALCULATED DATA/35X, 26HPERPENDICULAR DIPOLE(DE
1BYES),5X,28HAREA/MOLECULE(SQ, ANGSTROMS)/)
206  FORMAT (45X,F4,2,38X,F6,2)
500  CALL STOPPLOT

```

Security Classification		DOCUMENT CONTROL DATA - R & D	
(Security classification of title, body of abstract and indexing annotation must be entered when the overall report is classified)			
1. ORIGINATING ACTIVITY (Corporate author)		2a. REPORT SECURITY CLASSIFICATION	
Naval Research Laboratory Washington, D.C. 20390		Unclassified	
		2b. GROUP	
3. REPORT TITLE			
CDC 3800 COMPUTER PROGRAM TO CALCULATE AND PLOT THE COMPRESSIBILITY, MODULUS OF COMPRESSIBILITY, AND PERPENDICULAR COMPONENT OF THE MOLECULAR DIPOLE MOMENT OF FILMS SPREAD ON LIQUID SURFACES			
4. DESCRIPTIVE NOTES (Type of report and inclusive dates)			
Interim report; work is continuing on the problem.			
5. AUTHOR(S) (First name, middle initial, last name)			
Laura A. Halper			
6. REPORT DATE		7a. TOTAL NO. OF PAGES	7b. NO. OF REFS
December 31, 1969		7	
8a. CONTRACT OR GRANT NO.		9a. ORIGINATOR'S REPORT NUMBER(S)	
NRL Problem C02-10		NRL Report 7003	
b. PROJECT NO.			
RR 001-01-43-4751			
c.		9b. OTHER REPORT NO(S) (Any other numbers that may be assigned this report)	
d.			
10. DISTRIBUTION STATEMENT			
This document has been approved for public release and sale; its distribution is unlimited.			
11. SUPPLEMENTARY NOTES		12. SPONSORING MILITARY ACTIVITY	
		Department of the Navy (Office of Naval Research), Washington, D.C. 20360	
13. ABSTRACT			
<p>A program of interest to surface chemists has been developed for use on NRL's Control Data Corporation 3800 computer. The program reads film pressure, surface potential, and molecular area data. Tabulated values are returned for the compressibility, modulus of compressibility, and perpendicular component of the molecular dipole moment, along with corresponding molecular areas. Library routines for the Calcomp Plotter are used to plot the compressibility vs area/molecule, modulus of compressibility vs area/molecule, and perpendicular component of the molecular dipole moment vs area/molecule.</p>			

14. KEY WORDS	LINK A		LINK B		LINK C	
	ROLE	WT	ROLE	WT	ROLE	WT
Computers Surface properties Compressibility Modulus of compressibility Dipole moment Monomolecular films						